

NAME

pyFAI-recalib – Deprecated

DESCRIPTIONusage: pyFAI-recalib [options] **-i** ponifile **-w** 1 **-c** calibrant.D imagefile.edf

Calibrate the diffraction setup geometry based on Debye–Sherrer rings images with a priori knowledge of your setup (an input PONI–file). You will need to provide a calibrant or a "d–spacing" file containing the spacing of Miller plans in Angstrom (in decreasing order). Calibrants available: Ni, CrOx, NaCl, Si_SRM640e, Si_SRM640d, Si_SRM640a, Si_SRM640c, Si_SRM640b, Cr2O3, AgBh, Si_SRM640, CuO, PBBA, alpha_Al2O3, quartz, C14H30O, cristobaltite, Si, LaB6, CeO2, LaB6_SRM660a, LaB6_SRM660b, LaB6_SRM660c, TiO2, ZnO, Al, Au or search in the American Mineralogist database: <http://rruff.geo.arizona.edu/AMS/amcsd.php> The **—calibrant** option is mandatory !

positional arguments:

FILE List of files to calibrate

optional arguments:**-h, —help**

show this help message and exit

-V, —version

show program's version number and exit

-o FILE, **—out** FILE

Filename where processed image is saved

-v, —verbose

switch to debug/verbose mode

-c FILE, **—calibrant** FILE

Calibrant name or file containing d–spacing of the reference sample (MANDATORY, case sensitive !)

-w WAVELENGTH, **—wavelength** WAVELENGTH

wavelength of the X–Ray beam in Angstrom. Mandatory

-e ENERGY, **—energy** ENERGYenergy of the X–Ray beam in keV ($hc=12.398419292\text{keV}\cdot\text{\AA}$).**-P** POLARIZATION_FACTOR, **—polarization** POLARIZATION_FACTORpolarization factor, from **-1** (vertical) to **+1** (horizontal), default is None (no correction), synchrotrons are around 0.95**-i** FILE, **—poni** FILE

file containing the diffraction parameter (poni–file). MANDATORY for pyFAI–recalib!

-b BACKGROUND, **—background** BACKGROUND

Automatic background subtraction if no value are provided

-d DARK, **—dark** DARK

list of comma separated dark images to average and subtract

-f FLAT, **—flat** FLAT

list of comma separated flat images to average and divide

-s SPLINE, **—spline** SPLINE

spline file describing the detector distortion

-D DETECTOR_NAME, **—detector** DETECTOR_NAME

Detector name (instead of pixel size+spline)

-m MASK, **—mask** MASK

file containing the mask (for image reconstruction)

-n NPT, **--pt NPT**
file with datapoints saved. Default: basename.npt

--filter FILTER
select the filter, either mean(default), max or median

-l DISTANCE, **--distance DISTANCE**
sample-detector distance in millimeter. Default: 100mm

--dist DIST
sample-detector distance in meter. Default: 0.1m

--poni1 PONI1
poni1 coordinate in meter. Default: center of detector

--poni2 PONI2
poni2 coordinate in meter. Default: center of detector

--rot1 ROT1
rot1 in radians. default: 0

--rot2 ROT2
rot2 in radians. default: 0

--rot3 ROT3
rot3 in radians. default: 0

--fix-dist
fix the distance parameter

--free-dist
free the distance parameter. Default: Activated

--fix-poni1
fix the poni1 parameter

--free-poni1
free the poni1 parameter. Default: Activated

--fix-poni2
fix the poni2 parameter

--free-poni2
free the poni2 parameter. Default: Activated

--fix-rot1
fix the rot1 parameter

--free-rot1
free the rot1 parameter. Default: Activated

--fix-rot2
fix the rot2 parameter

--free-rot2
free the rot2 parameter. Default: Activated

--fix-rot3
fix the rot3 parameter

--free-rot3
free the rot3 parameter. Default: Activated

--fix-wavelength
fix the wavelength parameter. Default: Activated

--free-wavelength
free the wavelength parameter. Default: Deactivated

--tilt Allow initially detector tilt to be refined (rot1, rot2, rot3). Default: Activated

--no-tilt
Deactivated tilt refinement and set all rotation to 0

--saturation SATURATION
consider all pixel>max*(1-saturation) as saturated and reconstruct them, default: 0 (deactivated)

--weighted
weight fit by intensity, by default not.

--npt NPT_1D
Number of point in 1D integrated pattern, Default: 1024

--npt-azim NPT_2D_AZIM
Number of azimuthal sectors in 2D integrated images. Default: 360

--npt-rad NPT_2D_RAD
Number of radial bins in 2D integrated images. Default: 400

--unit UNIT
Valid units for radial range: 2th_deg, 2th_rad, q_{nm}^{-1} , q_A^{-1} , r_mm. Default: 2th_deg

--no-gui
force the program to run without a Graphical interface

--no-interactive
force the program to run and exit without prompting for refinements

-r MAX_RINGS, --ring MAX_RINGS
maximum number of rings to extract. Default: all accessible

-k, --keep
Keep existing control point and append new

The main difference with pyFAI-calib is the way control-point hence DebyeSherrer rings are extracted. While pyFAI-calib relies on the contiguity of a region of peaks called massif; pyFAI-recalib knows approximately the geometry and is able to select the region where the ring should be. From this region it selects automatically the various peaks; making pyFAI-recalib able to run without graphical interface and without human intervention (**--no-gui** and **--nointeractive** options). Note that 'pyFAI-recalib' program is obsolete as the same fonctionnality is available from within pyFAI-calib, using the 'recalib' command in the refinement process. Two option are available for recalib: the numbe of rings to extract (similar to the **-r** option of this program) and a new option which lets you choose between the original 'massif' algorithm and newer ones like 'blob' and 'watershed' detection.